

1 **ACPSEM position paper: Commissioning, and quality assurance of Magnetic Resonance**
2 **Imaging Linear Accelerators**

3
4 Luis Munoz¹, Reza Alinaghi Zadeh², Stephen Gibson³, Sandra Fisher², Nick Cook⁴,
5 Christopher Colyer¹, Zoe Moutrie^{5,6,7} on behalf of the ACPSEM Magnetic Resonance Imaging
6 Linac Working Group (MRILWG)

- 7
8 ¹GenesisCare Adelaide, SA, Australia
9 ² Austin Health, Melbourne, VIC, Australia
10 ³Townsville University Hospital, Townsville, QLD, Australia
11 ⁴Christchurch Hospital, Christchurch, New Zealand
12 ⁵ South Western Sydney Cancer Services, New South Wales Health, NSW Australia
13 ⁶ South West Sydney Clinical Campuses, University of New South Wales, Sydney, NSW,
14 Australia
15 ⁷ Ingham Institute for Applied Medical Research, Liverpool, NSW, Australia

16
17 Contents

18 **1. ABSTRACT2**
19 **2. INTRODUCTION2**
20 **3. SCOPE OF PRACTICE3**
21 **4. DEFINITIONS AND ABBREVIATIONS.....3**
22 **5. MACHINE OVERVIEW4**
23 **6. WG POSITION ON ACCEPTANCE AND COMMISSIONING OF AN MRI-**
24 **LINAC7**
25 **7. RECOMMENDATIONS ON QA DEVICES AND THEIR REQUIREMENTS.. 27**
26 **8. WORKING GROUP POSITION ON PERIODIC QA29**
27 **9. FUTURE DEVELOPMENTS AND CHALLENGES IN MRI-LINAC**
28 **TECHNOLOGY36**

29	10. SUMMARY	36
----	--------------------------	-----------

30	11. ACKNOWLEDGMENTS	36
----	----------------------------------	-----------

31

32 1. Abstract

33 The Magnetic Resonance Imaging Linac Working Group (MRILWG) present a position
34 statement on the commissioning and quality assurance (QA) tests for linear accelerators
35 coupled with Magnetic Resonance Imaging. The core objective of the MRI-Linac quality
36 assurance (QA) sub-group was to curate a set of critical performance tests to assist physicists
37 in establishing and maintaining a safe and effective treatment program.

38 The commitment to a vendor neutral approach was made to delineate recommendations
39 towards site preparation, commissioning assessments, QA tests and their frequency. The
40 foundational and longitudinal studies referenced in this endeavour emphasized a broad
41 spectrum of sources to provide a comprehensive guidance. Tests presented aim to reflect
42 clinical use patterns and intend to be sensitive and relevant to detecting errors related to the
43 specific use of an MRI-Linac. A certified ROMP is responsible for authorising return of the
44 radiation therapy equipment to clinical use following any repair, adjustment, upgrade or
45 modification to the equipment that affects patient safety.

46 2. Introduction

47 The MRI-Linac presents an advancement in the precise delivery of radiotherapy. Utilizing MR
48 imaging, there is enhanced distinction of soft tissue, with the added capability of procuring
49 dynamic and functional tissue data in real-time. An improved understanding of the tumour
50 microenvironment [1] provides the potential for true treatment adaptation, leading to further
51 escalation of dose to target tissues and optimized organ at risk and normal tissue sparing [2].

52 The integration of new technologies presents opportunities and challenges, spanning facility
53 planning, acceptance and commissioning and quality assurance (QA). To understand the
54 impact of an MR environment on beam generation and its interaction within and around a
55 patient, physicists determine which tests are pertinent at the time of acceptance and
56 commissioning and on-going QA. In this pursuit, references were made to foundational studies
57 like Roberts et al [3], longitudinal studies pertaining from both Elekta Unity (Elekta AB,
58 Stockholm, Sweden) and ViewRay MRIdian (MRIdian™, ViewRay Inc., Cleveland, OH,

59 USA), have been incorporated for this purpose [4-7]. Further insights from Woodings et. al.
60 [8] and the ESTRO-ACROP consensus opinion [9] augmented this work to ensure a
61 comprehensive guidance.

62 This position paper aims at:

- 63 - describing a comprehensive set of vendor agnostic recommendations on acceptance and
64 commissioning tests, supplemented by justification.
- 65 - detailing site preparation considerations.
- 66 - listing current QA devices
- 67 - recommendations on routine QA tests for daily, monthly and yearly frequencies.

68 As commercial MRI-Linac platforms evolve to add functionality and features, requirements
69 for acceptance, commissioning and QA should be reviewed in light of these changes.

70 3. Scope of practice

71

72 The purpose of this position statement is to address the commissioning and quality assurance
73 (QA) of MRI-Linacs. In view of the critical nature of acceptance, commissioning, and the
74 ongoing management of the QA program, specialized knowledge and training is required [10].

75 As such the working group recommend that the certified ROMP assumes the ultimate
76 responsibility for overseeing and executing the appropriate completion of tests, evaluations,
77 and assessments. They will also provide guidance on QA procedures and protocols, ensuring
78 that they adhere to best practices and regulatory standards [11]. Delegation of specific tests to
79 non-certified ROMPs should be completed with specific training, oversight and review of QA
80 results [12].

81 4. Definitions and Abbreviations

82

83	AAPM	American Association of Physicists in Medicine
84	ACPSEM	Australian College of Physical Scientists and Engineers in Medicine
85	ACR	American College of Radiology
86	ARPANSA	Australian Radiation Protection and Nuclear Safety Agency
87	DAT	Device Acceptance Test
88	DIMP	Diagnostic Imaging Medical Physicist
89	DSV	Diameter of Spherical Volume

90	EPID	Electronic Portal Imaging Device
91	ERE	Electron Return Effect
92	ESE	Electron Streaming Effect
93	FFF	Flattening Filter-Free
94	FOV	Field of View
95	IMRT	Intensity-Modulated Radiotherapy
96	MLC	Multi-leaf Collimator
97	MRIL	Magnetic Resonance Imaging Linac
98	NCRP	National Council on Radiation Protection & Measurements
99	PSQA	Patient Specific Quality Assurance
100	QA	Quality Assurance
101	RF	Radiofrequency
102	ROMP	Radiation Oncology Medical Physicist
103	TPS	Treatment Planning System

104 5. Machine overview

105 5.1 Current commercial machines

106 Currently two commercial systems are available to the Australian and New Zealand market,
107 the Elekta Unity (Elekta AB, Stockholm, Sweden) and the ViewRay MRIdian (ViewRay
108 Technologies Inc, Cleveland, OH, USA) Figure 1A and 1B respectively. Details of their
109 characteristics are summarised in Table 1. Both employ an S-band linac delivering flattening
110 filter free (FFF) photon beams perpendicular to the static magnetic field. Delivery is completed
111 via step-and-shoot intensity-modulated radiotherapy (IMRT) on both platforms that facilitate
112 adaptive treatment workflows.

113 The MR hardware and technology architecture at the heart of these devices overcome the
114 engineering challenge presented to adequately isolate both systems from each other and find
115 synergy to leverage their composite benefit. To achieve the challenges of high spatiotemporal
116 resolution for real time imaging that satisfies the requirements for adaptation [13]. This is
117 dependent on having good gradient performance, along with high slew rates that influence the
118 minimum attainable TR and TE for imaging [14].

119 Each MRI-Linac design approaches the integration of systems differently using a combination
120 of active and passive magnetic shielding. For the Elekta Unity, active shielding that isolates
121 the MRI system from the linac component is used to provide a low-field toroid for the linac
122 beam and sensitive components [15], and a central 15 cm gap in coils plus shimming is

123 employed to deliver the photon beam [16]. The MRIdian utilises split coils leaving a 28 cm
124 gap and linac components are mounted in ferromagnetic compartments around a ring of
125 ferromagnetic shields forming a magnetically shielded volume around the linac [17].

126 Another notable difference between the two systems is that the Unity employs the Agility
127 multileaf collimator (MLC) while the MRIdian uses a double stack and double focus MLC
128 without additional jaws, reducing the effective leaf width to half the physical width and further
129 reducing an interleaf leakage. Refer to Zhang et al [18] and Latifi et al [15] for further details
130 on each respective system. The two clinical systems provide translational treatment couches,
131 the Elekta Unity provides only longitudinal movements in the direction of the magnetic field,
132 whereas the ViewRay MRIdian system allows lateral, longitudinal and vertical translations,
133 allowing for corrections due to initial patient alignment.

134



135 Figure 1 – A) Elekta Unity and B) ViewRay MRIdian

136 5.2 First Generation Systems

137 5.2.1 Australian MRI-Linac

138 The Australian MRI-Linac is a prototype system that can deliver 4 and 6 MV photons through
139 a Varian Millennium MLC mounted on a rail system employing an open bore 1 T magnet [19].

140 The magnet is a split superconductive magnet with the ability to have the in-line or

141 perpendicular beam field orientation which is similar to the Aurora design [20].



142

143 Figure 2 - Australian MRI-Linac

144 5.2.2 Aurora System

145 The Aurora RT (MagnetTx Oncology Solutions, Canada) system [21] employs a 6MV linac
146 and 0.6 T MR, where the linac can be positioned between open MR planes (perpendicular), or
147 alternatively through a central opening of one of the planes (in-line) (Figure 3).



148

149 Figure 3 - Aurora MagnetTx system

150

151

152

153

154 Table 1. Configuration details of current MRI-Linac systems

155
156

Feature	Elekta Unity	MRIdian ViewRay	Australian MRI	Aurora MagnetTx
Nominal Energy (MV)	7 FFF	6 FFF	4 & 6 FFF	6 FFF
B0 strength (T)	1.5	0.35	1.0	0.55
SAD (cm)	143.5	90	190 - 330	122
MLC	Single stack 160 leaves	Double stack 138 leaves	Single stack 120 leaves	Single stack 120 leaves
MLC speed (mm/s)	86	40	-	35
Max field size (cm)	57.4 x 22 cm ²	27.4 x 24.1 cm ²	30 cm ² to 50cm ²	28.5 cm x 28.5 cm ²
Bore diameter (cm)	70	70	82	110W x 60H
Magnet type	Closed superconductor	Split superconductor	Open superconductor	Open Room temperature MR
Orientation	Perpendicular	Perpendicular	Inline and perpendicular	Inline
Delivery method	Step and shoot IMRT	Step and shoot IMRT	Step and shoot IMRT	Sliding window Step and shoot IMRT
360-degree delivery	No	Yes	Yes * with patient rotation	Yes
MRI Characteristics	gradient strength (GS) of 34 mT.m ⁻¹ and slew rate of 120 T.m ⁻¹ .s ⁻¹	GS 10 mT.m ⁻¹ and a slew rate of 200 T.m ⁻¹ .s ⁻¹	???	GS 45 mT.m ⁻¹ and a slew rate of 200 T.m ⁻¹ .s ⁻¹

157
158
159

160 **6. WG position on acceptance and commissioning of an MRI-Linac**
 161 MRI-Linac device acceptance tests are conducted with a vendor and customer component. Site
 162 physicists should participate during this process to ensure a level of accuracy and consistency
 163 commensurate with their local equipment. Recommended commissioning tests found in Table
 164 3 find synergy with ESTRO-ACROP recommendations detailed by Tanadini-Lang et al. [9],
 165 Woodings et al. (Unity) [8] and Valdenaire et. al. (MRIdian) [22]. Hybrid tests provided by
 166 Tijssen et al assess the interactions between linac and MRI system. Test descriptions are not
 167 provided in subsequent sections however, justification and expanded descriptions on some of
 168 the recommended tests are provided. Sections will be lettered accordingly.

169

170 Table 2 – Vendor agnostic recommended commissioning tests endorsed by the MRI linac
 171 working group.

172

Test/Check Category	Test/Check Description	Tolerances	References
A. Site preparation			
	Radiation Survey - Preliminary survey completed post installation and first beam on	As per local regulatory requirements	
	Final survey post beam tuning and machine calibration	As per local regulatory requirements	
MR survey	Complete magnetic fringe field assessment		[23]
Acoustic survey	Validation of acoustic insulation of RF cage/bunker and headphones/ hearing protection		[24]
MRI safety	Installation of metal detectors ferromagnetic detection system		[9],[24]
MRI infrastructure	Safe installation of quench pipe		
MRI influence on surrounding devices	Beam profile and output stability on adjacent linacs with MRI-linac gantry rotation pre and post magnet ramp	± 1% (TG142 tolerances)	[9, 22]
B. Acceptance tests (vendor performed)	Safety – Inhibit systems including emergency off, audio-visual, two-way communication system, interlocks, door, lights, beam stability	As per vendor specifications As per local regulatory requirements	
	Coordinate systems and data integrity	As per vendor specifications	[8, 25]
	Helium fill		
	Radiation isocentre – beam alignment and locus	± 1 mm	
	MV panel rigidity, alignment, pixel scale, isocentre and image quality	As per vendor specifications	[8]
	Beam quality	± 1%	[26]
	Dose output without gantry variation - MRL systems without specific TPS characterisation	≤ 1% (CCT usually 0.98 - 1.02)	
	Dose output with gantry rotation - MRL systems with specific TPS characterisation	± 2%	
	Dose rate stability with/without gantry rotation	± 2%	IEC 976 specify that the monitor chamber should have less than 2% variation to dose rate?

	MU linearity	2% > 5 MU	TG 142 and IEC 976 and IEC 977
	MU reproducibility	COV < 0.5	<IEC 976>
	Beam profile with gantry rotation	± 1%	[26]
	Beam limiting device calibration - MLC only or and MLC and jaw calibration	± 1 mm	[26]
	Gating and beam hold functionality	As per vendor specifications	
	Gating and beam hold latency	As per vendor specifications	
	MR to MV alignment	± 1 mm	TG284
	Individual coils and channels assessment	As per vendor specifications	
	Effect of MV beam on MR image quality	As per vendor specifications	
	Effect of gantry on MR image quality, and B0 homogeneity	As per vendor specifications	
	MR geometric accuracy	± 1 mm	
	System config checks and backups	As per vendor specifications	
C. Mechanical	Gantry angle calibration – rotation and readout	< 0.2°± 1°	
	Couch calibration, orthogonality and alignment	± 1 mm	TG 142 IEC 976
	MLC and jaw calibration, orthogonality, and sag	± 1 mm	IEC 976 AAPM TG 142
	QA support system calibration and alignment	± 1 mm	
	Alignment of all isocentres (laser, MV panel, beam, MRI)	± 1 mm	TG 142
D. Beam data collection		As per vendor requirements	
	Beam data required for creation of, and/or comparison against, beam model		
	Additional data collected for site specific TPS commissioning and ongoing routine QA.	As per clinical protocols and local regulatory requirements	
E. Dosimetry	Monitor chamber output constancy, stability, accuracy, and precision	± 1%	IEC 976
	MU Timer accuracy and system latency for gating	± 1% or ± 50 ms	IEC 976
	Beam quality and baseline	As per clinical protocols	TG 142

	Dose rate stability, linearity, reproducibility – short and long term	$\pm 1\%$	IEC 976
	Cryostat/ high density MR element characterisation	As per vendor specifications	
	Output with gantry angle	$\pm 1\%$	TG 142
	Flatness and Symmetry	$\pm 2\%$	TG 142
	Beam profile stability with gantry rotation	$\pm 1\%$	TG 142
	MLC and jaw transmission	$< 1\%$	TG119
	Couch transmission	As per clinical protocols	
	Coil transmission	As per vendor specifications	
	Immobilisation device commissioning	As per clinical protocols	
F. Ancillary Imaging	MV imager Central pixel location and subsequent geometric accuracy based on pixel size and mag factors	$\pm 1\%$	
	MV imager image quality	As per clinical protocols	TG 142
G. MRI Scanner	Cage and RF interference map	$\geq 100\text{dB}$ from Marlin 1.5T	
	MR to MV alignment and baselining, including multiple gantry angles	$\pm 1\text{ mm}$	
	MR uniformity and SNR measurement	acquired using all three coils (body, torso, H&N) in the transversal, sagittal, and coronal planes using the NEMA protocol for analysis	
	Magnetic field drift (B0 stability)	$< 1\text{ ppm/day}$ during $< 0.25\text{ ppm / day}$ for first 1-2 months operation	TG 284
	Transmitter and Gain Calibration	No visible artefacts Manual transmit gain within 5% of automatic Manual centre frequency within 10 Hz of automatic	TG 284
	Transmitter Gain Stability	Vendor specified minimum amplitude, frequency, and phase stability levels unless otherwise agreed upon	TG 284
	Magnetic field homogeneity (B0)	0.5 ppm volume root mean square (VRMS) across a 35cm DSV or as specified by MRI manufacturer across a specified DSV	TG 284

	Gradient Non-linearity	≤1 mm (within 10 cm radial distance of isocentre) ≤2mm (<20 cm radial distance away from magnet isocentre)	TG 284
	RF coil evaluation		TG 284
	External laser offset from MR isocentre	≤1 mm (where essential for patient positioning)	TG 284
	Table alignment with B0	$0 \pm 0.3^\circ$	TG 284
	Informatics/ connectivity/ Data transfer/ orientation	Site specific	TG 284 & TG 248
	MR spurious noise assessment	As per clinical protocols	[27]
	MR geometric distortion assessment and baselining	Acquire baseline for routine QA ≤ 2mm across 25 cm FOV	TG 284
	Effect of linac states on MR Image quality	As per vendor specifications	[27, 28]
	MR image quality and distortion with gantry rotation	As per vendor specifications	[29]
	MR image quality and distortion with MV beam on	As per vendor specifications	[29]
	MR image quality and distortion with MLC movement (single stack and double stack)	As per vendor specifications	[29]
	MR image quality and distortion with jaw (secondary MLC bank) movement	As per vendor specifications	[29]
	MR spatial integrity in cine mode	As per vendor specifications	[9]
	MR contrast and MR marker validation and protocol development	As per vendor specifications	
	Motion management assessment (gating)	As per vendor specifications	
	DWI and quantitative MRI assessment	As per vendor specifications	
H. TPS	Connectivity and acceptance testing to record and verify system	As per vendor specifications	
	Coordinate systems and data integrity	As per vendor specifications	[8]
	Basic field validation including output and symmetry	$\pm 2\%$	
	Heterogeneous field validation	As per clinical protocols	
	Couch and MR coil modelling validation	As per clinical protocols	
	Complex and simple clinical case validation	As per clinical protocols	
	Ancillary patient aid modelling e.g., headboard, masks, wing board	As per clinical protocols	
	ERE/ESE/ EF modelling in simple and anthropomorphic cases	As per clinical protocols	

	Density conversion method accuracy and equivalency during adaptations	As per clinical protocols	[9]
	DVH validation	As per clinical protocols	
	Density information layering (where appropriate) accuracy	As per clinical protocols	
	Contouring tool equivalency	As per clinical protocols	
	Validation and equivalency of adaption techniques	As per clinical protocols	
	End to end system tests, including dosimetry of simple and complex plans	As per clinical protocols	
	Clinical planning protocol or template commissioning	As per clinical protocols	
	Plan reconciling tool validation	As per clinical protocols	
I. End-To-End	End to end testing for clinical workflows	As per local regulatory requirements	
J. Patient plan specific QA	Secondary MU check program commissioning for simple and complex fields	As per clinical protocols	
	Measurement-based QA tool commissioning e.g., ArcCheck, phantom with film	As per clinical protocols	
Adaptive end-to-end test	Phantom studies for specific adaptive workflows	As per clinical protocols, suggested γ (5%/2mm), 10% threshold, >90% pixels passed	[8]
Non-adaptive end-to-end	Full end-to-end testing static targets	Townsville Paper – MRIL, Powers et al	
	Full end-to-end testing including motion management		
	Geometric fidelity and geometric consistency of all modes of acquisition of scanner i.e. Cine mode, navigated scan, localizer scan, fast and slow scans	4DMR sim/4DCT sim	R.A.
J. Audit	Absolute dosimetry audit by regulatory body or nearby MR linac centre	As per local regulatory requirements	

174 **6A MRI-Linac shielding and site preparation**

175 An MRI-Linac shares commonalities with conventional linear accelerators for radiation
 176 shielding with the added complexity of an MR environment. Optimal site setup should make
 177 provision to consider network requirements, RF shielding and isolation, quench vent and
 178 exhaust, excessive mechanical vibrations (whether steady state or transient) acoustic
 179 management and magnetic shielding which can impact an MRI-Linac system.

180 Site preparation guides are provided by vendors for standard configurations, the working group
 181 recommends close consultation with vendors during the planning process. Design aspects that
 182 improve safety and general site planning are detailed by Hu et. al. [30]. Of note is the need for
 183 dedicated storage for QA and immobilization equipment to limit misuse of unsafe equipment
 184 used for conventional linac treatment.

185 Internationally accepted shielding design protocols can be used to determine the most suitable
 186 material and attenuation level required with minor modifications. The MRI-Linac WG
 187 recommend NCRP 151 [31] as the values for tenth value layer (TVL) are more conservative.
 188 The differences as applicable to shielding between conventional and MRI-Linac are
 189 summarised in Table 3, each item for consideration is explained in detail in the paragraphs that
 190 follow. When considering neutron emission, standard photon shields within MRI linacs are
 191 adequate to safeguard against these neutrons at every available energy configuration.

192

193 *Table 3 2 Comparison of the pertinent machine characteristics for radiation shielding between*
 194 *conventional linacs and MR linac vendors.*

Section	Machine Characteristic	Parameter Affected	Conventional Linac	Elekta Unity	ViewRay MRIdian
T1	Source to isocentre distance	d_{SAD}	1000 mm	1435 mm	900 mm
	Primary collimator angle	<i>Primary barrier width</i>	27.8°	8°	18°
T2	Maximum field size	<i>Primary barrier width</i>	40 cm x 40 cm	57.4 cm x 22.0 cm	27.4 cm x 27.1 cm

	Maximum field area	F	1600 cm ²	1263 cm ²	743 cm ²
T3	Primary beam transmission	D_t	100%	0.5% at covers (2130 mm from iso)	Assume 100%
T4	Average leakage	L_f	0.1%	0.15%	0.1%
	Nominal dose rate	D_0	6-24 Gy/min	7 Gy/min	6 Gy/min

195

196 T1. Primary Barrier Transmission and Source to Isocentre Distance

197 In C-arm linacs, which typically have a distance of 1000 mm between the source and the
 198 isocentre, the source-to-isocentre distance is vitally important. Current MRI-Linac systems do
 199 not adhere to this standard distance. As all distances are normalized to 1 m in NCRP 151, the
 200 source-to-axis distance (dSAD) must be explicitly included. As an example, the calculation for
 201 primary barrier transmission would be as follows:

$$202 \quad B_{\text{pri}} = \frac{P \left(\frac{d_{\text{pri}}}{d_{\text{SAD}}} \right)^2}{WUT}$$

203

204 The primary beam transmission in MRI-Linacs is affected by the presence of the magnet,
 205 cryostat, and beam blockers. Vendors provide values for maximum transmission, which are
 206 measured at the machine cover or at a distance. Due to the inverse square law correction to the
 207 point, these values can be misleading about the overall amount of shielding provided. We
 208 recommend scaling back the maximum transmission value to isocentre, allowing it to be used
 209 in calculations of primary barrier thickness by scaling the workload or isocentric dose rate.

210 For example, considering the Elekta Unity:

211

$$212 \quad \text{Transmission} = 0.5\% \times \frac{(1435 + 2130)^2}{1435^2} = 3.09\%$$

213

214 In spite of the fact that this value is significantly higher than the average leakage of the unit, it
 215 is still common to find primary barriers that are thicker than adjacent secondary barriers. While

216 the beam spectra have already been altered by transmission through the machine, it is advisable
217 to use both TVL1 and TVLe for determining the primary barrier thickness, as the specific
218 alterations to the beam spectrum are not well-documented.

219

220 **T2 Maximum field size**

221 Conventional accelerators must consider the maximum field size to be defined by the jaws and
222 MLCs at a collimator rotation that produces the greatest lateral dimension. As Unity and
223 MRIdian have fixed collimators this consideration is not needed and the field size in the
224 superior and inferior direction can be simply projected onto the primary barrier. It is also
225 common to define the maximum field size by the primary collimator angle. The correct
226 maximum field size is therefore limited by the lower value of these two approaches.

227

228 **T3 Primary barrier width**

229 MRI-Linacs absorb much of the small-angle patient scatter due to their high energy. Scattering
230 angles are typically limited to $>25^\circ$, reducing penetration and scatter fraction. Consequently,
231 the requirement for the primary barrier to intercept at least the 20° scatter line is no longer
232 crucial, providing the opportunity to use laminated or composite barriers.

233 It is still recommended to add 30 cm either side of the largest beam projection to determine the
234 width of the primary barrier. This projection is maximized at the intersection of the ceiling
235 barrier and the wall barrier, as described in NCRP 151.

236 Steel can be used to shield and reinforce the bunker, but its use must be approved by the vendor
237 to ensure that the magnet can be adjusted to accommodate the additional steel. Non-ferrous
238 materials are required in and around the magnet, and each vendor offers specific guidance in
239 this regard.

240

241 **T4 Leakage radiation**

242 Some MR linear accelerator manufacturers report a higher percentage leakage than is observed
243 in conventional linacs. Thus, it is important to include this explicitly in the equation for leakage
244 transmission. We recommend this be denoted as the factor L_f , such that the leakage equation is
245 written as:

246

$$B_L = \frac{Pd_L^2}{L_f WT}$$

247

248 **Conduits, Magnetic fields, ventilation, air-conditioning and acoustic management**

249 Radiation therapy equipment typically requires conduits to facilitate cabling and dosimetry tool
250 installation through radiation shielding. These conduits are designed to minimize scattered
251 ionizing radiation outside the bunker [32]. Similarly, RF shielding allows for small
252 "waveguides" within the RF cage to accommodate conduits. Their design ensures that the RF
253 shield always remains effective, however, any conductive cabling that passes through a
254 waveguide may introduce RF artifacts into the MR images. As a result, it is crucial to minimize
255 the likelihood of stray RF signals interfering with the quality of MRI images. During
256 commissioning works or daily quality assurance, it may not always be possible to fully close
257 these waveguides. To determine the optimal arrangement, we recommend collaboration
258 between physicists, MR physicists, and the vendor.

259 The impact of an MRI-Linac on the surrounding environment and vice versa is crucial for a
260 successful build. It is important to monitor any major facility changes, such as construction in
261 adjacent areas or large ferromagnetic sources like elevators or MRI systems (above and below),
262 as they may affect the static field and distortions. MR image quality may be adversely affected
263 by transient as well as steady-state mechanical vibrations. In addition, stray magnetic fields
264 must be taken into account when high field MRI-Linacs are located near other linacs [33].
265 Particularly after static magnetic field ramp-up and any scheduled (or unscheduled) ramp-
266 down, ensure flatness and symmetry of gantry position on adjacent linacs.

267 Heat dissipation and management are crucial to MRI-Linac operation. Helium systems require
268 tight control of cooling, not just for cryostats but also for control cabinets. In order to meet the
269 vendor's operational specifications, physicists should consult cooling engineers when
270 designing MRI-Linac bunkers. Insufficient cooling for helium or high ambient temperatures
271 may result in extreme humidity in cryostat cabinet rooms, resulting in excessive condensation
272 in ancillary rooms, control consoles, and bunkers. It is recommended to establish robust cooling
273 chain monitoring overseen by local subject matter experts who can rectify issues outside of
274 regular hours. To address any disruptions to the cooling chain or ambient conditions, clear
275 monitoring protocols and management plans are essential. Monitoring environmental
276 conditions, such as temperature, humidity, water ingress, and linac performance, is useful for
277 pre-emptive maintenance and outage monitoring.

278

279 Managing acoustics is another consideration. Certain MRI procedures can produce moderate
280 to high levels of acoustic noise, affecting patient comfort and posing safety concerns for staff

281 and patients [34-36]. Through proper bunker construction, acoustic isolation can be achieved
282 to reduce MRI bunker noise [37]. Acoustic noise levels should be confirmed by physicists and
283 local standards should be followed [24].

284

285 **Summary recommendations**

286 - Vendor Consultation: Engage closely with vendors during the planning phase to
287 understand specific requirements, ensuring the optimal performance of MRI-Linac
288 systems.

289

290 - Shielding & Primary Barriers: Adopt the NCRP 151 protocol for conservative shielding
291 design. Pay keen attention to the primary barrier transmission, especially concerning
292 the differences in source-to-isocentre distances. For MRI-Linacs, prioritize the
293 consideration of small-angle patient scatter absorption and adjust primary barrier width
294 accordingly.

295

296 - Technical Specificities: Understand the nuances of maximum field sizes, especially
297 with MRI-Linacs' fixed collimators. Remain vigilant about potential leakage,
298 particularly focusing on the L_f factor that indicates higher leakage in MRI-Linacs.

299

300 - Setup & Environment: Minimize stray RF signals' interference with MR image quality,
301 necessitating close collaboration with MR physicists and vendors. Monitor the
302 surrounding environment of the MRI-Linac, including adjacent construction and major
303 facility changes.

304

305 - Temperature & Acoustic Control: Collaborate with cooling engineers to ensure
306 consistent temperature regulation vital for MRI-Linac operation. Emphasize acoustic
307 isolation in bunker construction to mitigate MRI-procedure noises, adhering to local
308 standards.

309

310 **6B MRI-Linac acceptance testing**

311 **6.B.1 Safety**

312 Standard safety tests include electrical, mechanical and dosimetric. Tests and tolerances
313 designed by AAPM TG-142 updated in TG-198 and acceptance tests within IEC 60976 [38]

314 and IEC 60601 [39] remain relevant as references to establish baseline guidelines for
315 acceptance tests.

316 A robust MRI safety program designed referencing ACR recommendations will ensure staff
317 knowledgeable of the risks inherent in working in an MR environment [40]. The MRI-Linac
318 WG endorses the recommendations presented by the MRI-Linac Safety WG [24].

319 **6.B.2 System configuration and connectivity**

320 Acceptance testing and commissioning for DICOM imaging and connectivity to the hospital
321 PACS, record and verify system (RVS) and other required nodes is a paramount consideration
322 for real-time adaptive workflows. Tijssen et al. provide a list of system connection and
323 configuration tests recommended at the time of acceptance and commissioning [25]

324

325 **6C Mechanical**

326 **6.C.1 Radiation isocentre**

327 Mechanical isocentre is known to influence radiation isocentre in C-arm linacs, as collimator
328 and couch do not rotate in current commercial MRI-Linac offerings, the concept of mechanical
329 isocentre does not strictly apply. The validation of co-incidence of radiation and imaging
330 isocentres has been verified to achieve a 1 mm tolerance by several groups across both
331 commercial platforms [8, 41, 42].

332 The size of the radiation iso-centre can be assessed using a star-shot image using film
333 sandwiched between copper plates as detailed by Roberts [3] or Palacios et al. [43], or
334 alternatively using the MV portal dosimeter if one is provided. The position and the size can
335 then be determined using a Winston-Lutz test, where any possible sag found in the beam
336 limiting device [38]. Powers et al provide an overview of the linac commissioning tests
337 performed on Elekta Unity along with their results, however with an emphasis on how to
338 perform some of the tests when dedicated, specialized equipment is not available [44].

339

340 **6.C.2 Considerations for the beam limiting device**

341 MLC and jaw alignment with gantry rotation and treatment beam is a stalwart component of
342 machine characterisation. Each MLC bank should be independently verified for the positioning
343 and transmission, for the MRIdian system this is particularly important where there exists two
344 sets of MLCs stacked in a tessellated configuration. Both systems have fixed collimator
345 positions and validation of any possible rotation should be a consideration as detailed by
346 Woodings [8]. Users are to familiarise themselves with the vendor's jaw and MLC calibration

347 processes, in using available detectors in the clinic independently verify alignment, interleaf
348 leakage and intra-leaf leakage over full MLC bank and extent of travel [18, 45]. Picket fence
349 tests using film or EPID can be completed for verification, Tsuneda et. al describes workflows
350 in overcoming limitations with the Elekta Unity [46].

351

352 **6D. Scanned and non-scanned data considerations**

353 Considerations for beam scanning presented by AAPM-TG106 hold relevance for acquisition
354 of baseline data for acceptance and commissioning [47]. For the acquisition of relative
355 dosimetry, the effect of the magnetic field on a scanning detector within a scanning water tank
356 may vary with depth, off axis position and field size, this effect should naturally be lower for
357 low field strength systems like the ViewRay.

358 Key points of consideration for scanning under the influence of the static magnetic field
359 include: profile offsets, changes in the effective points of measurement [48] and the influence
360 of detector orientation [49]. The introduced shift in the effective point of measurement of
361 ionization changes affects linac calibration and commissioning of the treatment planning
362 system. Specifically, the vertical shift impacts beam quality and tissue maximum ratio
363 determination, and the lateral shift affects ion chamber usage.

364

365 Baseline profile data can be acquired using film, planar array, or portal dosimetry if available.
366 According to Roberts, chamber response can be variable when using closely packed ion
367 chambers in an array due to differences in average density around the chambers. It is common
368 to use planar arrays for relative profile measurements during quality assurance, so cross-
369 calibration against water tank data would be necessary to account for device sensitivity and the
370 B-field effect on profile measurements [3].

371 An additional consideration is the presence of air gaps in phantoms which can impact dosimetry
372 measurements for the MRI-Linac, as the MRI's magnetic fields may influence dose distribution
373 and accuracy due to altered electron return effect in the air-filled regions, one solution is to fill
374 the gaps with water [50] [51].

375

376 **Summary Recommendations:**

- 377 - During validation of the beam data and acquisition of reference fields for constancy
378 measurement, the working group emphasise checking with vendor requirements on
379 which detectors are recommended for scanned and non-scanned beam data. This is to

380 ensure data acquired during the acceptance and commissioning period best emulates
381 the data during the modelling process.

- 382 - Select detectors that are known to provide accurate results in the presence of a magnetic
383 field, such as ionization chambers and diamond detectors. Avoid using shielded diodes,
384 as they can produce misleading dose profiles.
- 385 - A single detector type should be used to acquire scanned and non-scanned data.
- 386 - Corrections for lateral shift in the beam profile should be applied after OPF are applied.
- 387 - If different detectors are required, OPF should be measured at the point of peak
388 intensity.
- 389 - Consider the influence of the magnetic field on detector positioning: The magnetic field
390 may cause shifts in the effective point of measurement for detectors. To ensure accurate
391 measurements, use an on-board MV portal imaging system or another suitable method
392 to verify the reproducibility of detector positioning within the magnetic field.

393 **6E. Dosimetry - Absolute and Relative dose measurements**

394 The QA working group endorses the use of the MRI-Linac WG – dosimetry paper for reference
395 and recommendations when completing absolute dosimetry [52]. Relative dosimetry
396 measurements will experience the same dependencies attributed to measurement within an MR
397 environment.

398

399 **6.E.1 Dosimetric characterization of patient support, immobilization and ancillary** 400 **imaging equipment**

401 A variety of immobilization and accessory devices are used to ensure positional repeatability
402 during treatment, these devices must be correctly represented in the TPS for safe adaptive
403 workflows. The working group recommends that transmission measurements of patient support
404 hardware including the table and immobilization devices should be checked as part of the
405 acceptance and commissioning process [53].

406 Hu et al [20] provides useful recommendations when considering such devices. For the
407 MRIdian system this includes a fibreglass couch top which is moveable and the indexing and
408 lateral placement of patients on top may impact dose both at the surface and at depths. For the
409 Elekta Unity this includes the high-density couch support struts, it is recommended by the
410 vendor to avoid treating through, however there is no interlock to prevent this occurring offline
411 or online in the dedicated TPS.

412 The working group recommends that the receiver coils radiation beam transmission should also
413 be verified at the time of commissioning for each coil available. Liney and Raaijmakers

414 foresaw and evidenced the effect of RF receiver coil impacting dosimetry respectively [54]
415 [55], [56]. Powers et al [44] points out that no work to date had been published on anterior coil
416 attenuation characterization for the Elekta Unity, nor does it form part of the device acceptance
417 test; rather, a factory default structure and relative electron density (RED) is applied in the TPS.
418 The authors of this paper recognise the early work by Hoogcarspel et al. [57] which showed a
419 decrease in dose of up to 2.2% as a result of the coil, suggesting that modelling of the dosimetric
420 impact of the coil must be considered in planning.

421 At the time of this paper, Elekta provided a model in Monaco of the standard body coil to
422 predict the dosimetric impact of the device. The position and orientation of these coils may not
423 be guaranteed for all vendors. As such, extensive validation of the coil modelling should be
424 performed with care taken to understanding the impact of height above patient surface,
425 longitudinal positioning above the isocentre and any tilting that may occur fraction to fraction.
426 The Powers investigation addressed this for the Elekta Unity coil, finding only a maximum
427 difference of only 0.5% between measured and calculated (in Monaco) attenuation across the
428 lateral extent of the coil; however, the out of field dose due to ESE from the coil was found to
429 be significant [58].

430 Radiation induced currents (RIC) may impact RF coils leading to image artefacts as detailed
431 by [59, 60]. While Buckley et al [61], Hoogcarspel et al [57] and Burke et al [59] demonstrate
432 methods for assessing the RIC using Fast Field Echo sequences and the ACR MRI accreditation
433 phantom, as some of these artifacts may depend on the positioning, while others are static
434 influences and may not vary with clinical use. The working group recommends understanding
435 these influences which may lead to uncertainties during motion monitoring online.

436 **6F. Ancillary imaging systems**

437 The implementation of a robust QA program can rely on ancillary imaging systems provided
438 with a platform [44]. The working group recommends the use of AAPM physics practice
439 guideline 2b, for commissioning tests pursuant with desired outcomes [62].

440

441 **6G. MRI-Linac considerations on MRI**

442 There are several factors arising from the presence of an MR imaging system that can affect
443 the accuracy of treatment planning, these need to be investigated and quantified for appropriate
444 QA. The primary focus of MR performance for treatment planning is geometric fidelity. The
445 Unity and MRIdian manage the fringe field in quite different ways: the Unity low-field toroid
446 is designed to minimise the fringe field around sensitive components, and any potential beam

447 steering issues introduced from gantry or MV imager sag have been investigated to be within
448 recommended tolerances provided in TG142 [41]. The MRIdian shielded ring contains large,
449 unevenly distributed ferromagnetic components which are isolated from affecting the fringe
450 field and static field homogeneity. Ginn et al. assessed the MRIdian and provides results to
451 inform appropriate planning target volume (PTV) margins for 0.35 T MRI-guided radiotherapy
452 [63].

453 **6.G.1 Linac Gantry rotation and magnetic field homogeneity**

454 The ring design employed by both commercial offerings allows for the potential modification
455 of the fringe field of the imaging magnet as the linac rotates around the patient. This may cause
456 rotation-dependent changes in the imaging field homogeneity, linac performance, and beam
457 steering. For systems where the linac is fixed or where the magnet rotates with linac rotation,
458 software active shimming can provide a viable solution to disturbances in the B_0 field
459 homogeneity, however this requires active monitoring on the part of the physicist to ensure
460 appropriate compensation is made for each image. As an indication of system performance,
461 magnetic field drift tests in the initial 2 months of acceptance is recommended by TG 284 [23].

462 **6.G.2 Image quality measures**

463 Image quality of the MRI system in the MRI-Linac requires constant assessment. Image
464 uniformity, signal-to-noise ratio and spatial resolution are usually assessed by vendor-provided
465 phantoms during routine QA. It is also recommended to check the performance of individual
466 coils and channels based on the vendor guidance, including monitoring of long-term stability
467 [64]. Most MRI manufacturers provide semi-automated analysis phantoms and tools to meet
468 NEMA standards.

469 The ACR MRI quality control tests are well established and widely adopted in diagnostic MRI.
470 The same action levels and frequencies are recommended to be performed on the MR-Linac
471 system for sequence specific assessment [25].

472 ACR recommended parameters to assess include:

- 473 • high contrast resolution,
- 474 • slice thickness accuracy,
- 475 • slice position accuracy,
- 476 • image intensity uniformity,
- 477 • signal ghosting,
- 478 • low contrast detectability, and
- 479 • signal-to-noise ratio.

480 RF coil placement over treatment areas significantly impacts image quality. It is crucial to
481 understand the position and height of moveable RF coils above the patient, as signal
482 degradation increases with distance. Image quality metrics should be assessed for all height
483 variations, tilt, and longitudinal alignment over the intended target. Routine testing often
484 assumes fixed height, position, and minimal tilt. Systematic investigation of these factors, such
485 as assessing the impact of 2 cm shifts on image quality or the effect of key anatomy positioning
486 near the central imaging plate, is recommended. Lee et al describes a method to assess the
487 impact of coil tilt on image quality for reference [65].

488 **6.G.3 Geometric Distortion assessment**

489 No magnet is perfect, and the presence of a patient or phantom further disrupts the static field's
490 homogeneity. Precession frequency determines spatial encoding and is directly related to local
491 field strength. System-level geometric distortions in MRI result from static magnetic field
492 inhomogeneity, gradient magnetic field nonlinearity, and patient-level chemical shifts and
493 susceptibility distributions. System-level distortion increases with radial distance from the MRI
494 isocentre, while susceptibility distortion increases with magnetic field strength at interfaces
495 [67]. When commissioning the MRI-Linac system, both static magnetic field inhomogeneity
496 and gradient nonlinearity should be assessed individually and baselined by medical physicists.
497 They should also be checked during set frequencies and particularly following gradient
498 calibration and services. This measurement can be done using reverse readout of gradient
499 polarity technique, or by isolating the gradient nonlinearity distortion from B_0 homogeneity
500 [23].

501 In addition, increasing gradient strength (i.e., increasing readout frequency bandwidth) can
502 reduce susceptibility distortion and chemical shift effects. As recommended by AAPM TG284
503 [23], the combined geometric distortion in an MRI must not exceed 1 mm in a 20 cm diameter
504 spherical volume (DSV) and 2 mm in a 40 cm diameter spherical volume. Walker et al. [66]
505 propose a vendor neutral method for MRIGRT geometry distortion assessment. This allows
506 physicists to assess any distortion that may be present on MRI sequences between MRI
507 simulation and an MRI-Linac which could impinge on the effectiveness of the simulation
508 process. Ensure geometric fidelity testing is appropriate and meets the tolerances used in
509 treatments with off-axis targets or small organs.

510 Additional, sequence-dependent distortion is caused by induced eddy-currents resulting from
511 rapid switching of gradients and gradient nonlinearities, which increase with distance from

512 isocentre. Non-homogeneity and gradient nonlinearity can be measured and are, to some
513 extent, compensated for by shimming, gradient compensation, and software correction. These
514 compensations, whilst usually sufficient in diagnostic radiology may give residual distortions
515 that could severely impact the accuracy of RT [67].

516 **6.G.4 Sequence Assessment**

517 Since MRI images of the MR-Linac systems are the primary images for the adapted treatment
518 of the patient, it is highly recommended that all MRI sequences in the MR-Linac system be
519 evaluated for image quality and geometric distortion to provide better estimation for dosimetric
520 uncertainty resulting from the MRI system. Also, the orientation of MRI images for the
521 sequences transferred to the TPS must be checked during the commissioning process using a
522 phantom with directional differences to confirm correct orientation is retained [25].

523 **6.G.5 MRI to MV isocentre**

524 Similar to all modern linear accelerators with onboard imaging system, MRI-Linac systems
525 require characterization and minimisation of the offset between imaging isocentre and radiation
526 isocentre. Each MRL vendor provides dedicated phantoms and processes for the MR-to-MV
527 isocentre check. Baseline MR-to-MV isocentre information is acquired during commissioning
528 and is frequently checked as part of the routine QA process [3].

529 **6.G.7 Diffusion weighted and quantitative imaging**

530
531
532 Monitoring tumour response during a course of treatment and adaptively modifying the
533 treatment plan based on tumour biological feedback may represent a new paradigm for
534 radiotherapy [68]. Different parameters like longitudinal relaxation rate (T1), transverse
535 relaxation rate (T2), diffusion-weighted imaging (DWI) and apparent diffusion coefficient
536 (ADC) have the potential to provide clinical findings. Since the accuracy of the results for both
537 quantitative measures and ADC values can vary based on the magnet model, sequence
538 implementation and magnet characterization, [69, 70] it is recommended to baseline and QA
539 with recommended phantoms and sequences to isolate the system variation from patient
540 response. Furthermore, as the design of each of the currently available systems differs from
541 that of typical diagnostic systems, such as the split gradients in the 1.5T system and the low
542 field of the 0.35 T system, recommendations exist for DWI scanning parameters such as b-
543 values for both systems [71]. There are several studies on developing quality assurance process
544

545 for the quantitative measures on clinically available MR-Linacs [5, 68], QA of the technical
546 performance of MRI for quantitative imaging is recommended to ensure metrics and factors
547 affecting results are due to a physiological response and not measurement variability [72].

548

549 **6H. MRI-Linac WG position on Treatment Planning System**

550 **6.H.1 Modelling ancillary components**

551 For precise dosimetry during MRI-Linac commissioning, it's essential to characterize the
552 ancillary components dosimetrically [44]. Initial modelling helps establish clear planning
553 guidelines. Components like moveable imaging coils, headphones, and cabling that may
554 influence radiation exposure are already well modelled, as highlighted by Powers[44, 58].

555

556 MRI-Linac systems pose unique challenge in imaging immobilization and support systems
557 [30]. These include fibreglass or carbon fibre couches, bolus, compression belts, breast boards,
558 head and neck masks, foam pillows, and vacuum bags. There are difficulties in positioning
559 these support systems consistently for each fraction, as well as MRIs cannot show these
560 components. To circumvent this, some centres have adopted custom virtual models in their
561 TPS [73] or implemented MRI-visible marking systems for daily delineation. Considering their
562 dosimetric influences, alternative modelling methodologies, such as CTs or virtual CTs, are
563 pivotal. Physicists must be adept at understanding the implications of these systems across
564 diverse clinical sequences.

565 **6.H.2 Less familiar sources of radiation and accounting for them in the TPS**

566 The imaging magnet alters dose distribution in the patient and ancillary support devices,
567 potentially reducing dose delivered to the treatment volume or increasing dose to non-target
568 organs [74]. Readers are referred to [75, 76] and the MRWG Dosimetry paper [52] for more
569 information. The working group recommends testing the ability of a TPS to account for
570 electron return effect and electron streaming effect. It is important to manage these influences
571 carefully in the beam modelling and to understand them thoroughly for routine clinical
572 planning [74]. Hall et al demonstrates both effects in the same clinical setting using
573 perpendicular MRI-Linacs [77].

574

575 **6I. End to End commissioning**

576 Commercial vendors have provided several devices for end-to-end verification of online
577 adaptive therapy workflows, these are detailed in Table 4. Quasi-3D detectors like the Sun
578 Nuclear ArcCheck (Sun Nuclear Inc, Melbourne, FL, USA) or PTW Ruby (PTW, Freiburg,
579 Germany) can be used as an end-to-end validation tool, and the use of dynamic motion to test

580 adaptation workflows is also recommended [78]. Before making clinical decisions based on
581 device measurements, users should verify the device's accuracy [79], especially since
582 measurement systems can be affected by magnetic fields, as highlighted by Ellefson et al. [59].

583

584 When MR-conditional or MR-safe commercial equipment isn't available, using in-house, non-
585 commercial solutions is recommended [80, 81]. Without a quasi-3D device, Powers et al. detail
586 an end-to-end commissioning process for IMRT treatment on the Unity system using a PTW
587 1500MR 2D array, Gafchromic EBT3/EBT-XD film, and the EPID. They assessed individual
588 beam segments and analysed composite deliveries using film or detector arrays. For end-to-
589 end tests, they used a water phantom with 3D-printed elements mimicking different target sizes
590 [44].

591

592 True 3D dosimetry is possible with radiochromic polymer gel, providing solutions for machine
593 QA and end-to-end verification [82, 83]; the use of gel dosimetry for 4-D verification of
594 accumulated dose using gel dosimetry is similarly feasible [84, 85]. The advent of modular
595 phantoms for end-to-end verification of treatment systems is well-documented, however the
596 inclusion of gel dosimetry looks to be a promising validation tool for end-to-end [86] and
597 routine QA [87] with interesting results.

598 **6.I.1 Motion management assessment**

599 Intrafraction motion from respiratory, musculoskeletal, cardiac and gastrointestinal systems is
600 a known issue in modern radiotherapy and there are different imaging techniques and processes
601 recommended in guidelines and publications [88, 89]. It is recommended to use a motion
602 phantom with a known waveform, frequency and amplitude for testing the motion managed
603 imaging, Cine imaging, 4DMRI process and identify and characterize the system limitations
604 [23]. Some vendors already provide the necessary respiratory gated and breath-hold treatment
605 workflows to process tumour tracking. We also anticipate ECG gated workflows which are
606 currently in their development phase and will be available in future.

607

608 **6J. System audit**

609 The MRI-Linac QA working group recognise the added benefit of independent end-to-end
610 validation for treatment processes. Especially in emerging technologies, where no formal
611 auditing organisation offers end-to-end evaluation with the treatment processes available in
612 your clinic we suggest working with nearby or established MIREL clinics to provide onsite

613 independent measurement-based testing. Audit programs by ARPANSA and IROC are also in
614 development to ensure future recourse for level 1B and level III audits.

615

616 7. Recommendations on QA Devices and their requirements

617 7.1 Influences of MR on dosimetry tools

618 While it is not the purpose of this paper to discuss all aspects of dosimeter performance which
619 may be affected by the presence of a low or high strength magnet [75, 90-92], physicists should
620 familiarise themselves with the necessary corrections and considerations that should be made
621 before using each detector with an MRI-Linac. Hu et. al. provides four for consideration when
622 evaluating QA equipment for the MR environment: 1) projectile hazard effect due to
623 ferromagnetic components, 2) electronic components that can be damaged by the magnet or
624 interfered by the time-varying RF and gradient fields, 3) the impact on the measurement
625 accuracy by the magnet and 4) image artifacts and distortion caused by the device [30]. Roberts
626 et al provide guidelines for testing methods and detector limitations with perpendicular 1.5 T
627 MRI-Linac [3].

628 Table 4 lists QA devices used by consortium members for all kinds of QA. The working group
629 recommend users only use equipment that vendors establish as MR safe and are labelled
630 properly as defined by ASTM requirements [93].

631

632 7.2 MRI conditional equipment

633 The presence of a magnetic field requires the use of MR conditional equipment to perform
634 commissioning and quality assurance. Each site will have a unique assortment of QA
635 equipment types and models, including devices verified as MR Conditional by the supplier,
636 devices manufactured locally, and some unknown or legacy. In all cases, the safety of devices
637 should be assessed according to the criteria in published guidance by the ACPSEM MR-Linac
638 Safety Guidance [24]. It is recommended that an equipment register is established to record
639 MR Safety status and functionality checks.

640

641 7.3 Validation of software

642

643 In their article [94] Salomons et al recommend medical physicists in radiation oncology apply
644 strict quality control to their patient pathway software. Although vendors are careful to
645 communicate major changes to software and describe which bugs have been resolved, these
646 changes can (and do) have an impact on overall performance. We recommend keeping strict
647 records of the software versions used and reporting on them, as well as maintaining strong

648 communication with vendors about software updates. Additionally, and more importantly,
 649 having test environments where a pre-clinical release is tested and vetted assures a smooth
 650 transition between updates to software.

651 Table 3 - MRI equipment for routine quality assurance

QA Device Type	Device Manufacturer/Name	Potential Use cases
Planar array devices	SNC Daily QA MR SNC IC Profiler-MR PTW STARCHECK maxi ® MR PTW Octavius 1500 MR, 1600 MR	Output measurement Flatness and symmetry Field size
1D detectors	PTW Semiflex 3D MR PTW PinPoint® 3D MR PTW microdiamond PTW Semiflex Exradin A1SL Exradin A19 IBA cc13 IBA cc04 Farmer type NE2571, IBA FC65, PTW30013 OSLD	Absolute/relative dosimetry Beam data collection Output factors Beam model validation Dose calibration Surface dose
Water tank	PTW MP1 MR Manual Water Phantom PTW BEAMSCAN® MR In house solutions	Absolute dose measurements Relative dose measurements Beam data collection
Film	Ashland, EBT3, XD	Absolute dose measurements Relative dose measurements Beam data collection MLC/Jaw calibration
MRI QA devices	Modus MRID(3D) geometric distortion phantom Philips geometric distortion phantom PIQT phantom ACR MR-MV test phantom (Vendor Supplied) CIRS Abdo 4D phantom 3DONE phantom	Image quality Isocentre offset Motion assessment validation Sequence testing
Patient specific QA/Quasi 3D and true 3D	SNC ArcCHECK®-MR OCTAVIUS® 4D PTW Octavius 1500 MR, 1600 MR DELTA4 Gel dosimeters	Patient specific plan QA Output with gantry rotation
MV imager	Integrated MV imager	Image quality Phantom setup Isocentre MLC/jaw calibration
End to End testing	CIRS STEEV CIRS Freepoint CIRS IMRT thorax PTW Ruby CIRS Abdo 4D phantom CIRS ZEUS MRgRT phantom Elekta/Varian end to end phantoms including: <ul style="list-style-type: none"> • Gel • Inhouse designed phantoms • CIRS IMRT THORAX 	Patient specific QA Beam model validation Gating
Software	AQUA	Isocentre

	RIT PyLinac MU2Net RadCalc ClearCalc/ClearCheck Inhouse software Vendor Supplied software	MLC/jaw calibration Secondary dose check Data transfer validation
Others	MV alignment phantom LasVegas phantom Solid water Setup jigs MR safe Thermometer MR safe Rulers CT/MR compatible markers	MV isocentre MV images quality Absolute and relative dosimetry

652

653 **8. Working group position on Periodic QA**

654 **8.1 Routine QA**

655

656 Routine QA is designed to be a subset of commissioning tests that aim to verify any actionable
657 differences when comparing results to baseline [26]. The site physicist is challenged to ensure
658 comprehensive testing without redundancy, critical performance tests and their tolerances must
659 be relevant to clinical patterns and are likely to evolve as the MRI-linac platform continues to
660 mature. By implementing resources such as TG 100, the physicist can ensure a robust QA
661 program is established [95]. In Table 5-8, tests are grouped by class with recommendations for
662 optional execution. The certified Radiation Medical Oncology Physics Specialist (ROMP) is
663 responsible for implementing the appropriate tests.

664 Users are encouraged to use a testing cadence that builds confidence in the performance of the
665 system. The experimental techniques for the recommended QA tests will not be described at
666 length, where published work is recommended, the reader is encouraged to pursue a
667 comprehensive understanding of required workflows.

668 The importance of QA of adaptive workflows cannot be emphasised enough, Chen et al [96]
669 provide a comprehensive end-to-end for daily QA. This workflow ensures that performance
670 checks and communication are verified on the Unity system.

671 Table 5. Daily quality assurance with suggested optional tests

Category	Procedure	Tolerance	Optional	Reference
Dosimetry	X-ray output constancy (all energies)	3%		AAPM TG-142 [26] / TG-198
Mechanical	MLC performance	Visual inspection of picket fence.		
Mechanical	Laser alignment	2mm		
Safety	Door Interlock	Functional		
Safety	Warning lights	Functional		
Safety	MRI regional specific safety checks <ul style="list-style-type: none"> • Low oxygen sensor • Compressor chirp • Ferromagnetic detectors 	Functional		AAPM TG-284 [23]
Safety	<ul style="list-style-type: none"> • Patient duress alarm • AV system • Coils and patient accessories • MRI-Linac emergency trolley check 	Functional		
Imaging	<p>EPID Image quality - a SNR type test, especially if using EPID based PSQA or EPID for MLC tests</p> <p>MRI – National Electrical Manufacturers’ Association (NEMA) standards</p> <p>Signal to noise Scaling Transverse and coronal</p> <p>Uniformity Spatial linearity Slice Profile Spatial Resolution Central Frequency</p>		** ** ** ** **	AAPM TG-284 [23], ACR MRI QA [40, 97]
End to end	Routine patient for daily QA testing DICOM functionality and communication Adaptive online QA – Secondary MU check Testing workflow for scan, plan, transfer, QA and treat	Local tolerances in dedicated phantom		Chen et al. (Unity) [96]

672

673

674 Table 6. Weekly quality assurance with suggested optional tests

Category	Procedure	Tolerance	Optional	Reference
Dosimetry	Beam quality (TPR 20,10) with gantry angle	±1% (TG-142)	**	[44]
	Backup monitor chamber constancy	±2%	**	[98]
	Photon beam profile	±1% from baseline	**	
	MLC	Visual inspection of picket fence	**	
Imaging	MR to MV PIQT Geometric Distortion	Translations ±0.5 mm to baseline Rotations between the MR and MV coordinate systems: maximum rotation for each axis: ± 0.3 degrees Mean of the absolute value of the rotations about each axis: ≤0.2 degrees	** ** **	Roberts et al. [1]

675

676 Table 7. Monthly quality assurance with suggested optional tests.

Category	Procedure	Tolerance	Optional	Reference
Dosimetry	Output constancy	2%		AAPM TG-142 [26] TG 198
	Backup monitor chamber constancy	2%		
	Photon beam profile and energy measurement			
Mechanical	Setting vs radiation field for two patterns (non-IMRT)	2 mm		
Mechanical	Couch position accuracy	IMRT 2 mm SBRT/SRS 1 mm		
	Gantry angle accuracy and reproducibility	±1°		
	Picket fence/MLC position accuracy or leaf position accuracy	Visual observation of picket fence or 1 mm for IMRT		[26]

		field at 4 gantry angles		
	Radiation iso centre size WL	± 1 mm		[98]
Mechanical	Localization lasers	1mm		
Safety	Warning lights, interrupts and interlocks Safety inspection of console, magnet and equipment room Compressor chip Ferromagnetic detectors AV system MRI bore fan RF door Low oxygen sensor Helium level Receiver coil and accessories check Emergency power off switch	Functional		AAPM TG-284 [23]
Imaging	MRI			AAPM TG-284 [23], ACR MRI QA [40, 97]
	High contrast spatial resolution	≤ 1 mm		
	Low contrast detectability	Total number discernible spokes (for four slices) for fields < 3 T: 21 (0.3 T) to 36 (1.5 T), 40 for 3 T		
	RF coil testing	Functional		
	Large field of view 3D geometric distortion	Verify ≤ 2 mm across 25 cm FOV		
	Percent image uniformity	Head coil: $\geq 87.5\%$ for < 3 T		
	Signal to noise	Scan and B0 dependent		
	Spatial linearity	0.5%		
	Slice thickness accuracy	± 1.0 mm		
	Static field verification - Transverse B0 & B1 maps	ctr freq. < \pm <1ppm/day during acceptance,		

		<0.25ppm/day for first 1-2 months operation		
	Central frequency drift	Manufacturer specified		
	MR and MV coincidence/isocentre shifts	Translations: 1 mm to baseline		
	Scaling Transverse and coronal	Within 1 mm		
	Transmitter gain stability	± 5% from baseline		
	Cryostat check (Helium level check)	Against baseline		
	Table check	± 1 mm from isocentre		
	Hardware check	Functional and without damage		
	MV			
	Positioning/repositioning	≤ 1 mm		
	Low contrast visibility – Las Vegas phantom	Baseline		
	Image quality – uniformity, artifact check (MV)	Baseline		
End to End	End to End TPS QA, routine patient for Annual QA Second check MU software audit	Baseline		Chen et al. (Unity) [96]

677

678 Table 8. Yearly quality assurance with suggested optional tests

Category	Procedure	Tolerance	Optional	Reference
Dosimetry	Output constancy	1%		AAPM TG-142 [26] TG 198
	Backup monitor chamber constancy	2%		
	X-ray profile symmetry comparison from baseline	1%		
	Spot check of field size dependent output factors	2% for < 4x4cm ² 1% for ≥ 4x4cm ²		
	X-ray beam quality	1%		
	X-ray MU linearity (output constancy)	5% (2-4MU), 2% ≥ 5MU		
	X-ray output constancy vs dose rate	2%		
	X-ray output constancy vs gantry angle	1%		
Mechanical	Couch maximum travel range	1 mm		

		(0.3 T) to 36 (1.5 T), 40 for 3 T		
	RF coil testing	Functional		
	Large field of view 3D geometric distortion	Verify $\leq 2\text{mm}$ across 25 cm FOV		
	Percent image uniformity	Head coil: $\geq 87.5\%$ for $< 3\text{ T}$		
	Signal to noise	Scan and B0 dependent		
	Spatial linearity	0.5%		
	Slice thickness accuracy	$\pm 1.0\text{ mm}$		
	Static field verification - Transverse B0 & B1 maps	ctr freq. $< \pm < 1\text{ppm/day}$ during acceptance, $< 0.25\text{ppm/day}$ for first 1-2 months operation		
	Central frequency drift	Manufacturer specified		
	MR and MV coincidence/isocentre shifts	Translations: 1 mm to baseline		
	Scaling Transverse and coronal	Within 1 mm		
	Transmitter gain stability	$\pm 5\%$ from baseline		
	Cryostat check (Helium level check)	Against baseline		
	Table check	$\pm 1\text{ mm}$ from isocentre		
	Hardware check	Functional and without damage		
	MV			
	Positioning/repositioning	$\leq 1\text{ mm}$		
	Low contrast visibility – Las Vegas phantom	Baseline		
	Image quality – uniformity, artifact check (MV)	Baseline		
End to End	End to End TPS QA, routine patient for Annual QA Second check MU software audit	Baseline		Chen et al. (Unity) [96]

680 **9. Future developments and challenges in MRI-Linac** 681 **Technology**

682 Adaptive radiotherapy strategies in the light of new technologies poses a challenge to the
683 workforce. As detailed by Hogan et. al., training and credentialing requirements underpin the
684 safe and efficient delivery of treatment on the MRI-Linac [99]. As this technology continues
685 to evolve and mature, time to treat will invariably shorten. The use of automation and artificial
686 intelligence in plan preparation presents a viable option towards workflow optimization as
687 detailed by Künzel et al [100] and Spieler et al. [101].

688 There are several developments underway to the MRI-Linac platforms, including advanced
689 gated delivery [102], improved treatment workflow optimization, and helical delivery [103].
690 The working group concedes that as clinical platforms develop, pertinent considerations for
691 acceptance and commissioning specific to each commercial offering will evolve.

692 **10. Summary**

693 The MRI-Linac QA working group make the following summarised recommendations:

- 694 a. Execute an end-to-end system check whenever a new or updated procedure is
695 introduced.
- 696 b. Include patient stabilization and support device transmission as a part of the
697 commissioning process.
- 698 c. Adaptive radiotherapy workflows require specifically designed QA protocols to ensure
699 a comprehensive assessment is made.
- 700 d. Collaborative commissioning is encouraged, including all relevant craft groups to
701 ensure a complete capture of the clinical workflow is assessed comprehensively.

702 **11. Acknowledgments**

703 The quality assurance working group would like to express their sincere gratitude for the time
704 and effort in reviewing this position paper. All the insightful comments and suggestions were
705 invaluable in helping us refine and improve this resource. We would like to thank external
706 reviewers [ADD EXTERNAL REVIEWERS], Katrina Biggerstaff for their thorough and
707 thoughtful review. We appreciate the dedication and expertise that the MRI-Linac Working
708 Group have brought, and we are grateful for the opportunity to collaborate with such a
709 knowledgeable and committed team.

710 **References**

711

- 712 1. Datta, A., et al., *Delivering functional imaging on the MRI-linac: current challenges*
713 *and potential solutions*. Clinical Oncology, 2018. **30**(11): p. 702-710.

- 714 2. Hall, W.A., et al., *The transformation of radiation oncology using real-time magnetic*
715 *resonance guidance: A review*. European Journal of Cancer, 2019. **122**: p. 42-52.
- 716 3. Roberts, D.A., et al., *Machine QA for the Elekta Unity system: A Report from the Elekta*
717 *MR-linac consortium*. Medical physics, 2021. **48**(5): p. e67-e85.
- 718 4. Mittauer, K.E., et al., *Characterization and longitudinal assessment of daily quality*
719 *assurance for an MR-guided radiotherapy (MRgRT) linac*. Journal of applied clinical
720 medical physics, 2019. **20**(11): p. 27-36.
- 721 5. Subashi, E., A. Dresner, and N. Tyagi, *Longitudinal assessment of quality assurance*
722 *measurements in a 1.5 T MR-linac: Part II—Magnetic resonance imaging*. Journal of
723 Applied Clinical Medical Physics, 2022: p. e13586.
- 724 6. Subashi, E., et al., *Longitudinal assessment of quality assurance measurements in a 1.5*
725 *T MR-linac: Part I—Linear accelerator*. Journal of Applied Clinical Medical Physics,
726 2021. **22**(10): p. 190-201.
- 727 7. Wetscherek, A., et al., *Longitudinal Stability of MRI QA up to Two Years on Eight*
728 *Clinical 1.5 T MR-Linacs*. Frontiers in Physics, 2022: p. 526.
- 729 8. Woodings, S.J., et al., *Acceptance procedure for the linear accelerator component of*
730 *the 1.5 T MRI-linac*. Journal of Applied Clinical Medical Physics, 2021. **22**(8): p. 45-
731 59.
- 732 9. Tanadini-Lang, S., et al., *An ESTRO-ACROP guideline on quality assurance and*
733 *medical physics commissioning of online MRI guided radiotherapy systems based on a*
734 *consensus expert opinion*. Radiotherapy and Oncology, 2023: p. 109504.
- 735 10. Marsh, L., et al., *ACPSEM position paper on ROMP scope of practice and staffing*
736 *levels for magnetic resonance linear accelerators*. Physical and Engineering Sciences
737 in Medicine, 2023: p. 1-7.
- 738 11. Agency., A.R.P.a.N.S., *The Code for Radiation Protection in Medical Exposure (RPS*
739 *C-5)*, in *Radiation Protection Series (2019)*. 2019: Australia. p. 39.
- 740 12. Medicine, A.C.o.P.S.a.E.i., *The Role of Physicists, Scientists and Engineers in*
741 *Medicine in Australasia*. 2018, Australasian College of Physical Scientists and
742 Engineers in Medicine (ACPSEM) Advisory Forum: Australia.
- 743 13. Ipsen, S., et al., *Radiotherapy beyond cancer: target localization in real-time MRI and*
744 *treatment planning for cardiac radiosurgery*. Medical physics, 2014. **41**(12): p.
745 120702.
- 746 14. McRobbie, D.W., et al., *MRI from Picture to Proton*. 2017: Cambridge university press.
- 747 15. Kooreman, E.S., et al., *ADC measurements on the Unity MR-linac—A recommendation*
748 *on behalf of the Elekta Unity MR-linac consortium*. Radiotherapy and Oncology, 2020.
749 **153**: p. 106-113.
- 750 16. Raaymakers, B., et al., *Integrating a 1.5 T MRI scanner with a 6 MV accelerator: proof*
751 *of concept*. Physics in Medicine & Biology, 2009. **54**(12): p. N229.
- 752 17. Klüter, S., *Technical design and concept of a 0.35 T MR-Linac*. Clinical and
753 Translational Radiation Oncology, 2019. **18**: p. 98-101.
- 754 18. Zhang, K., et al., *Performance of a multileaf collimator system for a 1.5 T MR-linac*.
755 Medical Physics, 2021. **48**(2): p. 546-555.
- 756 19. Jelen, U., et al., *Dosimetric optimization and commissioning of a high field inline MRI-*
757 *linac*. Frontiers in oncology, 2020. **10**: p. 136.
- 758 20. Keall, P.J., M. Barton, and S. Crozier. *The Australian magnetic resonance imaging—*
759 *linac program*. in *Seminars in radiation oncology*. 2014. Elsevier.
- 760 21. Fallone, B.G. *The rotating biplanar linac—magnetic resonance imaging system*. in
761 *Seminars in radiation oncology*. 2014. Elsevier.
- 762 22. Valdenaire, S., et al., *Acceptance, commissioning and quality assurance of the*
763 *MRIdian®: Site experience and three years follow-up*. Cancer/Radiothérapie, 2023.

- 764 23. Glide-Hurst, C.K., et al., *Task group 284 report: magnetic resonance imaging*
765 *simulation in radiotherapy: considerations for clinical implementation, optimization,*
766 *and quality assurance.* Medical physics, 2021. **48**(7): p. e636-e670.
- 767 24. Cook, N., et al., *ACPSEM position paper: the safety of magnetic resonance imaging*
768 *linear accelerators.* Physical and Engineering Sciences in Medicine, 2023: p. 1-25.
- 769 25. Tijssen, R.H., et al., *MRI commissioning of 1.5 T MR-linac systems—a multi-*
770 *institutional study.* Radiotherapy and Oncology, 2019. **132**: p. 114-120.
- 771 26. Klein, E.E., et al., *Task Group 142 report: Quality assurance of medical accelerators.*
772 *Medical physics*, 2009. **36**(9Part1): p. 4197-4212.
- 773 27. Tijssen, R.H.N., et al., *MRI commissioning of 1.5T MR-linac systems – a multi-*
774 *institutional study.* Radiotherapy and Oncology, 2019. **132**: p. 114-120.
- 775 28. Liney, G.P., et al., *experimental results from a prototype high-field inline MRI-linac.*
776 *Medical physics*, 2016. **43**(9): p. 5188-5194.
- 777 29. Liney, G.P., et al., *Imaging performance of a dedicated radiation transparent RF coil*
778 *on a 1.0 Tesla inline MRI-linac.* Physics in Medicine & Biology, 2018. **63**(13): p.
779 135005.
- 780 30. Hu, Q., et al., *Practical safety considerations for integration of magnetic resonance*
781 *imaging in radiation therapy.* Practical radiation oncology, 2020. **10**(6): p. 443-453.
- 782 31. d'Errico, F., *Structural shielding design and evaluation for megavoltage x-and gamma-*
783 *ray radiotherapy facilities: NCRP Report No. 151 Published by: National Council on*
784 *Radiation Protection and Measurements, 7910 Woodmont Avenue, Suite 400/Bethesda,*
785 *MD 20814-3095, USA ISBN: 0-929600-87-8, 246 pp (2005) 100.00(Hardback); 80.00*
786 *(electronic pdf file, available at <http://ncrppublications.org>).* Radiation Protection
787 Dosimetry, 2006. **121**(3): p. 342-343.
- 788 32. Sung, J., et al., *Compact bunker shielding assessment for 1.5 T MR-Linac.* Scientific
789 Reports, 2022. **12**(1): p. 1-9.
- 790 33. Perik, T., J. Kaas, and F. Wittkämper, *The impact of a 1.5 T MRI linac fringe field on*
791 *neighbouring linear accelerators.* Physics and Imaging in Radiation Oncology, 2017.
792 **4**: p. 12-16.
- 793 34. Price, D.L., et al., *Investigation of acoustic noise on 15 MRI scanners from 0.2 T to 3*
794 *T.* Journal of Magnetic Resonance Imaging: An Official Journal of the International
795 Society for Magnetic Resonance in Medicine, 2001. **13**(2): p. 288-293.
- 796 35. Carr, C.M., et al., *Evaluation of hearing loss in young adults after exposure to 3.0 T*
797 *MRI with standard hearing protection.* The Journal of the Acoustical Society of
798 America, 2022. **151**(3): p. 1913-1921.
- 799 36. Glans, A., et al., *Health effects related to exposure of static magnetic fields and acoustic*
800 *noise—comparison between MR and CT radiographers.* European Radiology, 2022: p.
801 1-14.
- 802 37. Evans, J.B., P. JEACOUSTICS, and N. SOLUTIONS, *SOUND ISOLATION DESIGN*
803 *FOR A MAGNETIC RESONANCE IMAGING SYSTEM (MRI).*
- 804 38. Commission, I.E., *Medical Electrical Equipment: Medical Electron Accelerators-*
805 *Functional Performance Characteristics.* 1989: International Electrotechnical
806 Commission.
- 807 39. Commission, I.E., *Medical electrical equipment—Part 2.1. Particular requirements for*
808 *the basic safety and essential performance of electron accelerators in the range 1 MeV*
809 *to 50 MeV.* International Electrotechnical Commission, 2009. **3**: p. 132.3. 4.
- 810 40. Safety, E.P.o.M., et al., *ACR guidance document on MR safe practices: 2013.* Journal
811 of Magnetic Resonance Imaging, 2013. **37**(3): p. 501-530.
- 812 41. Hunt, J.R., et al., *Variation in isocentre location of an Elekta Unity MR-linac through*
813 *full gantry rotation.* Physics in Medicine & Biology, 2022. **67**(1): p. 015005.

- 814 42. Wen, N., et al., *Evaluation of a magnetic resonance guided linear accelerator for*
815 *stereotactic radiosurgery treatment*. Radiotherapy and oncology, 2018. **127**(3): p. 460-
816 466.
- 817 43. Palacios, M.A. and C. Brink, *QA of MR-linac*, in *Advances in Magnetic Resonance*
818 *Technology and Applications*. 2022, Elsevier. p. 169-191.
- 819 44. Powers, M., et al., *Commissioning measurements on an Elekta Unity MR-Linac*.
820 *Physical and Engineering Sciences in Medicine*, 2022: p. 1-17.
- 821 45. Latifi, K., R. Lotey, and V. Feygelman, *On the MLC leaves alignment in the direction*
822 *orthogonal to movement*. Journal of Applied Clinical Medical Physics, 2021. **22**(6): p.
823 268-273.
- 824 46. Tsuneda, M., et al., *Elekta Unity MR-linac commissioning: mechanical and dosimetry*
825 *tests*. 2023, Oxford University Press.
- 826 47. Das, I.J., et al., *Accelerator beam data commissioning equipment and procedures:*
827 *report of the TG-106 of the Therapy Physics Committee of the AAPM*. Medical physics,
828 2008. **35**(9): p. 4186-4215.
- 829 48. O'Brien, D.J., et al., *Relative dosimetry with an MR-linac: response of ion chambers,*
830 *diamond, and diode detectors for off-axis, depth dose, and output factor measurements*.
831 *Medical physics*, 2018. **45**(2): p. 884-897.
- 832 49. Pojtinger, S., et al., *Ionization chamber correction factors for MR-linacs*. Physics in
833 *Medicine & Biology*, 2018. **63**(11): p. 11NT03.
- 834 50. Woodings, S.J., et al., *Beam characterisation of the 1.5 T MRI-linac*. Physics in
835 *Medicine & Biology*, 2018. **63**(8): p. 085015.
- 836 51. Hackett, S., et al., *Consequences of air around an ionization chamber: Are existing*
837 *solid phantoms suitable for reference dosimetry on an MR-linac?* Medical physics,
838 2016. **43**(7): p. 3961-3968.
- 839 52. Begg, J., et al., *ACPSEM position paper: dosimetry for magnetic resonance imaging*
840 *linear accelerators*. Physical and Engineering Sciences in Medicine, 2023: p. 1-17.
- 841 53. Olch, A.J., et al., *Dosimetric effects caused by couch tops and immobilization devices:*
842 *report of AAPM Task Group 176*. Medical physics, 2014. **41**(6Part1): p. 061501.
- 843 54. Liney, G.P., et al., *MRI-linear accelerator radiotherapy systems*. Clinical Oncology,
844 2018. **30**(11): p. 686-691.
- 845 55. Ghila, A., B. Fallone, and S. Rathee, *Influence of standard RF coil materials on surface*
846 *and buildup dose from a 6 MV photon beam in magnetic field*. Medical Physics, 2016.
847 **43**(11): p. 5808-5816.
- 848 56. Raaymakers, B.W., A.J. Raaijmakers, and J.J. Lagendijk, *Feasibility of MRI guided*
849 *proton therapy: magnetic field dose effects*. Physics in Medicine & biology, 2008.
850 **53**(20): p. 5615.
- 851 57. Hoogcarspel, S.J., et al., *Characterization of the first RF coil dedicated to 1.5 T MR*
852 *guided radiotherapy*. Physics in Medicine & Biology, 2018. **63**(2): p. 025014.
- 853 58. Powers, M. and J. Baines, *Electron Streaming Effect Associated With the Elekta Unity*
854 *Anterior Imaging Coil*. Frontiers in Physics, 2022: p. 293.
- 855 59. Burke, B., et al., *Radiation induced current in the RF coils of integrated linac-MR*
856 *systems: The effect of buildup and magnetic field*. Medical physics, 2012. **39**(8): p.
857 5004-5014.
- 858 60. Zijlema, S.E., et al., *Design and feasibility of a flexible, on-body, high impedance coil*
859 *receive array for a 1.5 T MR-linac*. Physics in Medicine & Biology, 2019. **64**(18): p.
860 185004.
- 861 61. Buckley, J.G., B. Dong, and G.P. Liney, *Imaging performance of a high-field in-line*
862 *magnetic resonance imaging linear accelerator with a patient rotation system for fixed-*

- 863 gantry radiotherapy. *Physics and Imaging in Radiation Oncology*, 2020. **16**: p. 130-
864 133.
- 865 62. McCullough, S.P., et al., *AAPM MEDICAL PHYSICS PRACTICE GUIDELINE 2. b.:
866 Commissioning and quality assurance of X-ray-based image-guided radiotherapy
867 systems*. *Journal of Applied Clinical Medical Physics*, 2021. **22**(9): p. 73-81.
- 868 63. Ginn, J.S., et al., *Characterization of spatial distortion in a 0.35 T MRI-guided
869 radiotherapy system*. *Physics in Medicine & Biology*, 2017. **62**(11): p. 4525.
- 870 64. Tijssen, R.H.N., et al., *MRI commissioning of 1.5T MR-linac systems - a multi-
871 institutional study*. *Radiother Oncol*, 2019. **132**: p. 114-120.
- 872 65. Lee, H.B., et al., *Image quality assessments according to the angle of tilt of a flex tilt
873 coil supporting device: An ACR phantom study*. *Journal of Applied Clinical Medical
874 Physics*, 2021. **22**(5): p. 110-116.
- 875 66. Walker, A., et al., *Development of a vendor neutral MRI distortion quality assurance
876 workflow*. *Journal of Applied Clinical Medical Physics*, 2022: p. e13735.
- 877 67. Wang, D., et al., *Geometric distortion in clinical MRI systems: Part I: evaluation using
878 a 3D phantom*. *Magnetic resonance imaging*, 2004. **22**(9): p. 1211-1221.
- 879 68. Gao, Y., et al., *Distortion-free diffusion MRI using an MRI-guided Tri-Cobalt 60
880 radiotherapy system: sequence verification and preliminary clinical experience*.
881 *Medical physics*, 2017. **44**(10): p. 5357-5366.
- 882 69. Fedeli, L., et al., *Dependence of apparent diffusion coefficient measurement on
883 diffusion gradient direction and spatial position—A quality assurance intercomparison
884 study of forty-four scanners for quantitative diffusion-weighted imaging*. *Physica
885 Medica*, 2018. **55**: p. 135-141.
- 886 70. van Houdt, P.J., et al., *Phantom-based quality assurance for multicenter quantitative
887 MRI in locally advanced cervical cancer*. *Radiotherapy and Oncology*, 2020. **153**: p.
888 114-121.
- 889 71. Nardini, M., et al., *Tuning the optimal diffusion-weighted MRI parameters on a 0.35-T
890 MR-Linac for clinical implementation: A phantom study*. *Frontiers in Oncology*, 2022.
891 **12**: p. 6622.
- 892 72. Carr, M.E., et al., *Conformance of a 3T radiotherapy MRI scanner to the QIBA
893 Diffusion Profile*. *Medical Physics*, 2022.
- 894 73. Dunkerley, D.A., et al., *Clinical Implementational and Site-Specific Workflows for a
895 1.5 T MR-Linac*. *Journal of Clinical Medicine*, 2022. **11**(6): p. 1662.
- 896 74. Baines, J., M. Powers, and G. Newman, *Sources of out-of-field dose in MRgRT: an
897 inter-comparison of measured and Monaco treatment planning system doses for the
898 Elekta Unity MR-linac*. *Physical and Engineering Sciences in Medicine*, 2021. **44**(4):
899 p. 1049-1059.
- 900 75. Jelen, U. and J. Begg. *Dosimetry needs for MRI-linacs*. in *Journal of Physics:
901 Conference Series*. 2019. IOP Publishing.
- 902 76. de Pooter, J., et al., *Reference dosimetry in MRI-linacs: evaluation of available
903 protocols and data to establish a code of practice*. *Physics in Medicine & Biology*,
904 2021. **66**(5): p. 05TR02.
- 905 77. Hall, W.A., et al., *Magnetic resonance linear accelerator technology and adaptive
906 radiation therapy: An overview for clinicians*. CA: a cancer journal for clinicians, 2022.
907 **72**(1): p. 34-56.
- 908 78. Liu, X., et al., *End-to-end verification of an MR-Linac using a dynamic motion
909 phantom*. *Medical Physics*, 2021. **48**(9): p. 5479-5489.
- 910 79. Ellefson, S.T., et al., *An analysis of the Arc CHECK-MR diode array's performance for
911 ViewRay quality assurance*. *Journal of Applied Clinical Medical Physics*, 2017. **18**(4):
912 p. 161-171.

- 913 80. Kim, K., et al., *Low-cost MR-compatible pneumatic respiratory organ motion*
914 *simulator for development of MR-guided thermal therapy*. Medical Physics.
- 915 81. Munoz, L., et al., *An open source solution for an in-house built dynamic platform for*
916 *the validation of stereotactic ablative body radiotherapy for VMAT and IMRT*.
917 Australasian physical & engineering sciences in medicine, 2016. **39**(4): p. 957-964.
- 918 82. Baldock, C., et al., *Polymer gel dosimetry*. Physics in Medicine & Biology, 2010. **55**(5):
919 p. R1.
- 920 83. Elter, A., et al. *RSC: Gel dosimetry as a tool for clinical implementation of image-*
921 *guided radiotherapy*. in *Journal of Physics: Conference Series*. 2022. IOP Publishing.
- 922 84. De Deene, Y. and M. Wheatley. *Real time 4D Radiation Gel Dosimetry on the*
923 *Australian MRI-Linac*. in *Journal of Physics: Conference Series*. 2022. IOP Publishing.
- 924 85. De Deene, Y., et al., *Towards real-time 4D radiation dosimetry on an MRI-Linac*.
925 Physics in Medicine & Biology, 2020. **65**(22): p. 225031.
- 926 86. Nierer, L., et al., *Evaluation of an anthropomorphic ion chamber and 3D gel dosimetry*
927 *head phantom at a 0.35 T MR-linac using separate 1.5 T MR-scanners for gel readout*.
928 *Zeitschrift für Medizinische Physik*, 2022.
- 929 87. Kim, J.H., et al., *3D star shot analysis using MAGAT gel dosimeter for integrated*
930 *imaging and radiation isocenter verification of MR-Linac system*. Journal of Applied
931 Clinical Medical Physics, 2022: p. e13615.
- 932 88. Keall, P.J., et al., *The management of respiratory motion in radiation oncology report*
933 *of AAPM Task Group 76 a*. Medical physics, 2006. **33**(10): p. 3874-3900.
- 934 89. Li, H., et al., *AAPM Task Group Report 290: Respiratory motion management for*
935 *particle therapy*. Medical physics, 2022. **49**(4): p. e50-e81.
- 936 90. Begg, J., et al., *Experimental characterisation of the magnetic field correction factor,*
937 *for Roos chambers in a parallel MRI-linac*. Physics in Medicine & Biology, 2022.
938 **67**(9): p. 095017.
- 939 91. Madden, L., et al., *In-line MRI-LINAC depth dose measurements using an in-house*
940 *plastic scintillation dosimeter*. Biomedical Physics & Engineering Express, 2021. **7**(2):
941 p. 025012.
- 942 92. Lim, S.B., et al., *An evaluation of the use of EBT-XD film for SRS/SBRT commissioning*
943 *of a 1.5 Tesla MR-Linac system*. Physica Medica, 2022. **96**: p. 9-17.
- 944 93. Woods, T.O., *Standards for medical devices in MRI: present and future*. Journal of
945 Magnetic Resonance Imaging: An Official Journal of the International Society for
946 Magnetic Resonance in Medicine, 2007. **26**(5): p. 1186-1189.
- 947 94. Salomons, G. and D. Kelly, *Software safety in radiation therapy*. Journal of Medical
948 Physics/Association of Medical Physicists of India, 2013. **38**(1): p. 1.
- 949 95. Huq, M.S., et al., *The report of Task Group 100 of the AAPM: Application of risk*
950 *analysis methods to radiation therapy quality management*. Medical physics, 2016.
951 **43**(7): p. 4209-4262.
- 952 96. Chen, X., et al., *A daily end-to-end quality assurance workflow for MR-guided online*
953 *adaptive radiation therapy on MR-Linac*. Journal of Applied Clinical Medical Physics,
954 2020. **21**(1): p. 205-212.
- 955 97. Safety., A.C.o.M., et al., *ACR guidance document on MR safe practices: Updates and*
956 *critical information 2019*. Journal of Magnetic Resonance Imaging, 2020. **51**(2): p.
957 331-338.
- 958 98. Hanley, J., et al., *AAPM Task Group 198 Report: An implementation guide for TG 142*
959 *quality assurance of medical accelerators*. Medical physics, 2021. **48**(10): p. e830-
960 e885.

- 961 99. Hogan, L., et al., *Old dogs, new tricks: MR-Linac training and credentialing of*
962 *radiation oncologists, radiation therapists and medical physicists*. Journal of Medical
963 Radiation Sciences, 2022.
- 964 100. Künzel, L.A., et al., *First experience of autonomous, un-supervised treatment planning*
965 *integrated in adaptive MR-guided radiotherapy and delivered to a patient with prostate*
966 *cancer*. Radiotherapy and oncology, 2021. **159**: p. 197-201.
- 967 101. Spieler, B., et al., *Automatic segmentation of abdominal anatomy by artificial*
968 *intelligence (AI) in adaptive radiotherapy of pancreatic cancer*. International journal
969 of radiation oncology, biology, physics, 2019. **105**(1): p. E130-E131.
- 970 102. Crijns, S., et al., *Towards MRI-guided linear accelerator control: gating on an MRI*
971 *accelerator*. Physics in Medicine & Biology, 2011. **56**(15): p. 4815.
- 972 103. Kontaxis, C., et al., *Proof-of-concept delivery of intensity modulated arc therapy on the*
973 *Elekta Unity 1.5 T MR-linac*. Physics in Medicine & Biology, 2021. **66**(4): p. 04LT01.
974
975